

## **PRODUCT DATA SHEET**

Product Data Sheet:	Sheep Monoclonal Antibodies (SMAs) to	
	benzodiazepines (benzo)	
SMA Covered	benzo.1D7; 1G3, 1G5, 1H5, 5A7, 5D3, 5D11, 10G10	
Introduction and	The measurement of <b>benzos</b> is useful as part of general	
<b>Production Description</b>	drug screening. The benzo family of drugs is sizeable	
-	and human samples contain the parental drug and a	
	complex mix of metabolites.	

## **Manufacturing Specification (both SMAs)**

Parameter	Specification	Test Method
Buffer and pH	PBS +0.1% Sodium azide at 7.2-7.4	pH Meter
Antibody Concentration	1-3 mg/mL	A280 (Ec = 1.4)
Purity and Aggregates	Single Main Peak ≥90%	Size Exclusion HPLC
Immunological Activity	±20% biological activity of frozen reference standard	In house ELISA

## **IMMUNOGENS:**

Chlorazepate-3-CO-X	1G5, 1H5, 5D3, 5D11, 10G10
Lorazepam-3-CO-X	1G3
Temazepam-3-CO-X	1D7

## **ELISA DATA: Free Benzo and Metabolite Affinities**

Please e-mail info@bioventix.com for detailed ELISA data

Datasheet last revised Nov 2015