


25(OH) Vitamin D2/D3

(Mouse monoclonal Antibody anti Vitamin D2 / D3)



PRODUCT ANALYSIS SHEET

Cat 53.197.06

Clone:	LMBP7013CB
Batch:	B9H2P1
Size	2 mg / 352 µl
Format:	Purified , Liquid (Phosphate Buffer)
Concentration:	5.67 mg/ml in PBS 0.1M pH7.4
Preservative:	NaN3 (0.05%)
Purification:	Purified by Protein A affinity chromatography
Isotype:	Ig G1 kappa
Immunogen:	Vitamin D derivate – BSA
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specificity :	Recognizes natural Vitamin D2 and Vitamin D3
Application:	Elisa
Storage & Stability:	2-8°C/10 years
Update:	

Date & Signature

DIAsource ImmunoAssays S.A.

Tel +32 10 84 99 11 | Fax +32 10 84 99 90 | www.diasource.be | Technical Inquiries: tech.support@diasource.be
To order: Belgium Free Phone 0800 159 59 | BeLux Tel +32 (0) 10 84 99 00 | Fax +32 (0) 10 84 99 96 | ordering@diasource.be
France Free Phone 0800908 443 | Free Fax 0800 902 588 | ordering@diasource.be